

NoIR LaserShields®

For Medical Applications



#38W
Universal Fitover



#34 - Wrap w/ Rx-Insert



#35 - Wrap around

LOOK
SMART™

Hair Removal:

- 810nm: DI8
- 980/1064nm: VG3
- 755nm: LIA
- 755/1064nm: CYN

Skin/Vein/Tattoo:

- 532/1064nm: DBY
- 585-600nm: DYE
- 10,600nm: EC2
- 2940nm: ERB

Lipo:

- 980/1064nm: VG3
- 1320nm: IRD

Combinations:

- 532/755/1064nm: TRI
- 810/980nm: ML1

IPL:

- 190-1200nm: 2PL
- 190-1200nm: 3PL

Patient:

- IShield (reusable):
- SmartShield (disposables):

Short List - Call for full listing of Aesthetic, Surgical, Ophthalmic and Dental LaserShields.

NoIR LaserShields provide industry-defining protection and visibility with an eye on comfort and wearability. With more filters and more frame options for medical applications, LaserShields are the point where safety and comfort converge.

All LaserShields are made in the USA with a lifetime warranty.

NoIR Laser Co, LLC

6155 Pontiac Trail • South Lyon, MI • 48178
800.321.9744 • 734.769.5565 • 734.769.1768 fax

www.lasereyeprotection.com

